

## III PULMONARY SUPPLIES

### PREVENT II™ FILTER

Part Number 536713-002, 536713-005 (Pediatric), 536713-004 (FOT)  
(100 each per box)

The preVent II™ filter uses high quality electrostat® material that traps bacteria and viruses, preventing contamination to the patient, technician and equipment. With a unique oval mouthpiece design and lip grip, the filter is comfortable to use. Designed to fit directly on the preVent flow sensor, the smart design is leak-free yet can be rotated for proper patient alignment.

- Eliminates the need for an additional coupler
- Filter design provides low resistance and minimal deadspace
- Eliminates leaks for accurate results
- Meets or exceeds ATS/ERS standards and specifications
- Provides cost effective barrier filtration for infection control
- Made in the U.S.A.

### KITS

Kit 1: Part Number 536719-001 (100 per box)  
contains the preVent II™ filter and the disposable, padded nose clip.

Kit 2: Part Number 536723-001 (100 per box)  
contains the preVent II™ filter, the disposable padded nose clip and the disposable mouthpiece.

Kit 3: Part Number 536719-004 (100 per box)  
contains the FOT filter and the disposable, padded nose clip.

Kit 4: Part Number 536723-004 (100 per box)  
contains the FOT filter, the disposable padded nose clip and the disposable mouthpiece.


These recommended kits contain parts in a sealed bag which create an entire patient testing package. These kits ensure efficiency, maximize clinical use and promote patient satisfaction.



PEDIATRIC



FOT




PREVENT II

### PREVENT II® FILTER SPECIFICATIONS

- Bacterial Efficiency: >99.98%,
- Resistance: 0.1 cm H<sub>2</sub>O/L/sec @ 0.5 L/sec
- Viral Filtration Efficiency: >99.92%
- Materials:
  - K resin case
  - Electrostat® filtrate medium
- Dead space: 53 mL
- Recyclable, single patient use



KIT 1 &  
KIT 3



KIT 2 &  
KIT 4

