



N555BL
Blue

N555AQ
Aqua

N555VL
Violet

N555BP
Bright Pink

N555PK
Rose

N555WH
White

N555SL
Silver

N555BK
Black

BRIGHTNESS COMPARISON

Product	# Lights	White	Red	Red/White	Batteries
NeoGlo [®]	4 LEDs	65,360 lx	X	60,210 lx	1
	2 LEDs	36,790 lx	14,760 lx	X	
Veinlite [®] EMS PRO	8 LEDs	39,180 lx	X	X	2
	16 LEDs*	X	43,840 lx	X	
Wee-Sight	3 LEDs	X	14,610 lx	X	2

*Veinlite's red light setting consists of 12 orange LEDs and 4 red LEDs.
Measurement done with an Extech Instruments Light Meter with sensor model T020-001.