





N555*



N5536RD



N5536WH

BRIGHTNESS COMPARISON

Product	# Lights	White	Red	Red/White	Batteries
NeoGlo [®] N555*	4 LEDs	65,360 lx	X	60,210 lx	1
	2 LEDs	36,790 lx	14,760 lx	X	
NeoGlo Disp N5536RD	2 LEDs	X	16,816 lx	X	0 (USB)
NeoGlo Disp N5536WH	2 LEDs	29,093 lx	X	X	0 (USB)
Wee-Sight	3 LEDs	X	14,610 lx	X	2

Measurement done with an Extech Instruments Light Meter with sensor model T020-001.

*N555 available in eight colors. See sell sheet M631.

©2021 Neotech Products LLC. All rights reserved.