

FreeArm Muscle with Peg

For the GE Giraffe OmniBed


The FreeArm Muscle with Peg standardizes gravity syringe feeds and venting inside the incubator and warmer, then transitions to open cribs, chairs and home with the interchangeable clamp.

Benefits:


- Holds feeding syringe and adjusts feed rate with the bendable arm.
- Holds open syringe out of the way for venting the stomach, improving patient safety.
- Frees nurses hands, improving nursing workflow.




1. Loosen flower screw.



2. Insert arm into peg or clamp and tighten flower screw.



3. Insert peg into Giraffe OmniBed corner or attach clamp to a flat surface or rounded bar.



Product Specifications:

Item #	Description	Qty/Case (each)
FAMPEG001X	Feed + Vent Syringe Holder (X=G (green) or W (white))	1 single / 24 per case
CLIP001001	Single Patient Use Clip	1 single / 25 per case

To order or for more information, email info@freearmcare.com.