



Bendable TarryBoards™

Flexible Support, Gentle Stability

Key Features:

- Provides comfort for the patient while supporting the I.V. site
- Handy adjustment tab at either end of board to help bend
- Green end threads make identification easy versus other TarryBoards
- Covered in soft spun fabric
- Single Patient Use
- Not made with natural latex rubber (latex-free)
- Made in U.S.A

Bendable TarryBoards™ offer customizable limb support with a bendable aluminum core that holds its shape for secure IV site stabilization, and oversized foam to enhance comfort. Soft spun fabric ensures comfort, while green end threads and easy-adjust tabs make positioning simple and clear. Optional straps eliminate the need for tape directly on the skin.



White Bendable TarryBoards™ by Size:



Preemie

T 72640 B100	Bendable TarryBoard, Preemie (1" x 2 1/2"), Single Use	100/box
T 72640 TVS B100	Strapped Bendable TarryBoard, Preemie (1" x 2 1/2"), Single Use	100/box

Neonatal

T 72641 B100	Bendable TarryBoard, Neonatal (1 1/2" x 4"), Single Use	100/box
T 72641 TVS B100	Strapped Bendable TarryBoard, Neonatal (1 1/2" x 4"), Single Use	100/box

Pediatric

T 72642 B100	Bendable TarryBoard, Pediatric (2" x 5"), Single Use	100/box
T 72642 VS B100	Strapped Bendable TarryBoard, Pediatric (2" x 5"), Single Use	100/box

Child

T 72645 B100	Bendable TarryBoard, Child (3" x 7"), Single Use	100/box
T 72645 VS B100	Strapped Bendable TarryBoard, Child (3" x 7"), Single Use	100/box

Pediatric

T 72693 B100	Bendable TarryBoard, Child Plus (3" x 9"), Single Use	100/box
T 72693 VS B100	Strapped Bendable TarryBoard, Child Plus (3" x 9"), Single Use	100/box

Pink and Blue Bendable TarryBoards™ by Size:

Preemie

T 72740 B100	Bendable TarryBoard, Preemie (1" x 2 1/2"), Asst Pink and Blue, Single Use	100/box
T 72740 TVS B100	Strapped Bendable TarryBoard, Preemie (1" x 2 1/2"), Single Use	100/box

Neonatal

T 72741 B100	Bendable TarryBoard, Neonatal (1 1/2" x 4"), Asst Pink and Blue, Single Use	100/box
T 72741 TVS B100	Strapped Bendable TarryBoard, Neonatal (1 1/2" x 4"), Single Use	100/box

Pediatric

T 72742 B100	Bendable TarryBoard, Pediatric (2" x 5"), Asst Pink and Blue, Single Use	100/box
T 72742 VS B100	Strapped Bendable TarryBoard, Pediatric (2" x 5"), Asst Pink and Blue, Single Use	100/box

Additional ordering information available on www.dandlelionmedical.com