



# Dandle® ROO Lite

## Position and Support Babies During Phototherapy

### The Dandle ROO Lite:

- Is available in 4 sizes to support babies up to 3600 g.
- Provides the caregiver with “another pair of hands” to support and comfort the infant during procedures
- Makes an excellent disposable positioning support during messy procedures, eye exams, transport, or after admission

The Dandle ROO Lite is a cost-effective solution for providing babies with developmentally supportive positioning that mimics the womb.

The soft, stretchy, lightweight cotton fabric also allows phototherapy to penetrate through to the baby’s skin. The ROO Lite’s patented design enables full customization of fit and function to meet every baby’s unique needs.

The wings can be used one at a time, two at a time, or in combination with the pouch. Using the pouch alone provides calming containment from feet to shoulders, while allowing up to 74% of phototherapy light to penetrate to the skin. To accommodate a bili blanket inside the positioner, the ROO Lite can be used in an alternate orientation with the baby’s feet nestled into the head roll (shown). The ROO Lite works well with a variety of available supportive accessories, including adjustable brims to shield the eyes from light, gel pillows to comfort and cushion, and a multipurpose roll to support midline flexion.



### Ordering Information

D 11031 B12 | Dandle ROO Lite - XS <1000g. | Qty: 12/box  
D 11032 B12 | Dandle ROO Lite - S 1000–1800 g. | Qty: 12/box  
D 11033 B12 | Dandle ROO Lite - M 1800–2500 g. | Qty: 12/box  
D 11034 B12 | Dandle ROO Lite - L 2500–3600 g. | Qty: 12/box  
See website for part numbers for additional accessories..

U.S. Patents: 9,332,791, 9,474,669, 9,820,907

