

Thicker pediatric gel pillows shown with white cotton covers



Standard gel pillows shown with washable printed cotton covers



Dandle•LION® Gel Pillows

With Disposable, White or Printed Cotton Covers

Size	Description
<i>(Standard thickness: ~5/8")</i>	
3" x 6"	Small Pillow
6" x 6"	Square Pillow
6" x 9"	Also fits XS Dandle ROO Lite or ROO2
6" x 9"	Also fits S Dandle ROO Lite or ROO2
6" x 12"	Also fits M Dandle ROO Lite or ROO2
<i>(Pediatric or Adult thickness: ~7/8")</i>	
6" x 7"	Pediatric Gel Pillow
6" x 10"	Pediatric Gel Pillow
<i>(Pediatric or Adult thickness: ~1")</i>	
8" x 13"	Pediatric Gel Pillow

Dandle•LION soft gel-filled pillows are designed to provide a comfortable cushion to meet the support needs of micropreemie, preemie, and full term infants.

Gel pillows can help reduce the need for frequent repositioning and make immobilized babies more comfortable. They can also be used as a mattress for very low birthweight babies or as a prone roll for bigger babies.

- Three covers available: white cotton, printed cotton, and disposable non-woven
- Non-toxic, water based food gel
- Guaranteed for six months
- DEHP and BPA free; not made with natural latex rubber
- Single patient use (unless hospital validates their own cleaning process)



Dandle·LION® Gel Pillow Ordering Information



3" x 6"		
D 14558 GP B12	Gel Pillow w/ printed cotton covers - 3" x 6"	12/box
D 14558 GPC B25	Gel Pillow Replacement printed cotton covers - fits 3" x 6" pillow	25/box
D 12558 GPC B25	Gel Pillow Disposable white covers - fits 3" x 6" pillow	25/box
D 0958 GP B12	Gel Pillow w/ white cotton covers - 3" x 6"	12/box
D 0958 GPC B25	Gel Pillow Replacement white cotton covers - fits 3" x 6" pillow	25/box
6" x 6"		
D 14559 GP B12	Gel Pillow w/ printed cotton covers - 6" x 6"	12/box
D 14559 GPC B25	Gel Pillow Replacement printed cotton covers - fits 6" x 6" pillow	25/box
D 12559 GPC B25	Gel Pillow Disposable white covers - fits 6" x 6" pillow	25/box
D 0959 GP B12	Gel Pillow w/ white cotton covers - 6" x 6"	12/box
D 0959 GPC B25	Gel Pillow Replacement white cotton covers - fits 6" x 6" pillow	25/box
6" x 9"		
D 14510 GP B12	Gel Pillow w/ printed cotton covers - 6" x 9"	12/box
D 14510 GPC B25	Gel Pillow Replacement printed cotton covers - fits 6" x 9" pillow	25/box
D 12551 GPC B25	Gel Pillow Disposable white covers - fits 6" x 9" pillow	25/box
D 0951 GP B12	Gel Pillow w/ white cotton covers - 6" x 9"	12/box
D 0951 GPC B25	Gel Pillow Replacement white cotton covers - fits 6" x 9" pillow	25/box
6" x 11"		
D 14520 GP B12	Gel Pillow w/ printed cotton covers - 6" x 11"	12/box
D 14520 GPC B25	Gel Pillow Replacement printed cotton covers - fits 6" x 11" pillow	25/box
D 12552 GPC B25	Gel Pillow Disposable white covers - fits 6" x 11" pillow	25/box
D 0952 GP B12	Gel Pillow w/ white cotton covers - 6" x 11"	12/box
D 0952 GPC B25	Gel Pillow Replacement white cotton covers - fits 6" x 11" pillow	25/box
6" x 12"		
D 14530 GP B12	Gel Pillow w/ printed cotton covers - 6" x 12"	12/box
D 14530 GPC B25	Gel Pillow Replacement printed cotton covers - fits 6" x 12" pillow	25/box
D 12553 GPC B25	Gel Pillow Disposable white covers - fits 6" x 12" pillow	25/box
D 0953 GP B12	Gel Pillow w/ white cotton covers - 6" x 12"	12/box
D 0953 GPC B25	Gel Pillow Replacement white cotton covers - fits 6" x 12" pillow	25/box
6" x 7"		
D 14555 GP B12	Pediatric Gel Pillow w/ printed cotton covers - 6" x 7"	12/box
D 14555 GPC B25	Pediatric Gel Pillow Replacement printed cotton covers - fits 6" x 7" pillow	25/box
D 12555 GPC B25	Pediatric Gel Pillow Disposable white covers - fits 6" x 7" pillow	25/box
D 0955 GP B12	Pediatric Gel Pillow w/ white cotton covers - 6" x 7"	12/box
D 0955 GPC B25	Pediatric Gel Pillow Replacement white cotton covers - fits 6" x 7" pillow	25/box
6" x 10"		
D 14556 GP B12	Pediatric Gel Pillow w/ printed cotton covers - 6" x 10"	12/box
D 14556 GPC B25	Pediatric Gel Pillow Replacement printed cotton covers - fits 6" x 10" pillow	25/box
D 12556 GPC B25	Pediatric Gel Pillow Disposable white covers - fits 6" x 10" pillow	25/box
D 0956 GP B12	Pediatric Gel Pillow w/ white cotton covers - 6" x 10"	12/box
D 0956 GPC B25	Pediatric Gel Pillow Replacement white cotton covers - fits 6" x 10" pillow	25/box
8" x 13"		
D 19054 B6	Lg. Pediatric gel pillow w/ white cotton covers - 8" x 13"	6/box
D 19054 GPC B25	Lg. Pediatric gel pillow Replacement white cotton covers - fits 8" x 13" pillow	25/box
D 12554 GPC B25	Lg. Pediatric gel pillow Disposable white covers - fits 8" x 13" pillow	25/box