



TarryGel™ Pillows

With White Cotton Covers

Key Features:

- DEHP-free, BPA-free polyethylene bag made without natural latex rubber
- Water-based food-grade gel that is resistant to bacterial and fungal growth and guaranteed mold-resistant
- Pillow covers are made of a flame-retardant polyester blend material that can be washed and reused
- Single Patient Use
- Made in U.S.A

TarryGel pillows are designed to provide a comfortable cushion to meet the support needs of micropremie, preemie, and full-term infants.

Gel pillows can help reduce the need for frequent repositioning and make immobilized babies more comfortable. They can also be used as a mattress for very low birthweight babies or as a prone roll for bigger babies.

Ordering Information

T 92332 B12	Small Gel Pillow, 3" x 6", 100g fill	12/case
T 92335 B12	Medium Gel Pillow, 6" x 6", 235g fill	12/case
T 92338 B12	Large Gel Pillow, 6" x 9", 400g fill	12/case
T 92332 C B50	Small Gel Pillow Covers, 3" x 6", Fits all 3" x 6" pillows	50/case
T 92335 C B50	Medium Gel Pillow Covers, 6" x 6", Fits all 6" x 6" pillows	50/case
T 92338 C B50	Large Gel Pillow Covers, 6" x 9", Fits all 6" x 9" pillows	50/case

