



# Vinyl TarryBoards™ & Sleeves

Durable Design. Comfortable Support.

## Key Features:

- Provides comfort for the patient while supporting the I.V. site.
- Handy adjustment tab at each end
- Lined with a thick foam pad to conform gently to the limb
- Covered in DEHP-free vinyl fabric
- Single Patient Use
- Not made with natural latex rubber (latex-free)
- Made in U.S.A

Vinyl TarryBoards™ combine foam padding and a DEHP-free vinyl cover for lasting support and easy cleaning. With gentle limb conformity, handy adjustment tabs, and the same trusted construction as our other TarryBoards™, they're a dependable solution.

## Ordering Information

T 40736 B100	Vinyl TarryBoard, Neonatal (1 1/2" x 4"), Single Patient Use	100/box
T 40737 B100	Vinyl TarryBoard, Pediatric (2" x 5"), Single Patient Use	100/box
T 40738 B100	Vinyl TarryBoard, Child (3" x 7"), Single Patient Use	100/box
T 40739 B100	Vinyl TarryBoard Adult (4" x 9"), Single Patient Use	100/box
T 40742 B100	Vinyl TarryBoard Adult Large (4" x 14"), Single Patient Use	100/box
T 40740 B100	Vinyl TarryBoard Adult X-Large (4" x 17"), Single Patient Use	100/box
T 40741 B100	Vinyl TarryBoard Adult XX-Large (5" x 18"), Single Patient Use	100/box

